

CHEMICALS AND TOXIC CABINETS





CHEMICALS AND TOXIC CABINETS FOR STORAGE

In the Chemicals & Toxic Cabinets, the space is optimally utilized with removable wire baskets, pull-out trays, removable shelves or our DIN FlexModul®, which can be used where particularly high flexibility and a logistics system are required for storage.

The interior design system DIN FlexModul® is the result of intensive development work to achieve maximum flexibility, systematic and safe storage with good clarity and order.

Our modular tray / basket system can be used wherever optimized space and tidiness are required, and the system can be moved between cabinets, industrial racks and mobile transport units.

Cabinets are provided with standard locks or system locks that fit the desired locking plan.

The cabinets are designed to withstand intensive use.



⚠ HOLD RIDEN LUKET
KEEP SASH CLOSED



Storage for DIN modules

FlexModul® stands for modular flexibility and high quality for absolutely clean storage units, for 60 cm wide cabinet modules.

Separator compartment is easily attached to the ClipOn fastening module and can be adjusted in intervals of 10 and 20 mm.

Likewise, sliding room divisions can be offered in all positions.

The trays and baskets can of course be equipped with end stops for fixed withdrawals. It is also possible to add its own touch on the slopes, and baskets with text or company logo.


Advantage

- Separate rooms can be adapted to the content
- Integrated drawer stop
- Leakage secured
- Low weight
- ABS plastic, scratch resistant and shockproof
- Temperature resistant from -40 ° C to + 70 ° C
- Large selection of sharing options.







Chemical cabinets



Chemical cabinets - double



Chemical Cabin - Front Pullout





Fronts

- Highpressure laminate
- Compact laminate
- Linolium
- Corian




Glass cabinet fronts

- Natural anodised 20 mm aluminum frame with 110° hinges
- White powder coated 20 mm aluminum frame with 110° hinges
- Laminate 60 mm frame with 270° hinges
- Laminate 60 mm frame concealed 110° hinges



Handles

- Standard logistics aluminum handles with inclined surface for for example Dymo strips. Width 190 mm
- Handle in brushed stainless steel. Width 160 mm
- Handle in matt aluminium. Width 160 mm



Countertops

Material options

- Laminate
- Compact laminate
- Stainless steel
- Corian
- Glass

Laminate colors



Fronts

Materials

- Laminate
- Compact laminate
- Corian

Laminate colors



Plinths

Two standard colors




Laminate colors





Tabletop profiles

Front edge



Rear edge





Postform




Compact laminate with Alu upturned edge





Trespa + laminate with Alu upturned edge





Stainless steel, with bead



Epoxy Resin with bead



Polypropylen with bead



Ceramic with bead



Trays and baskets

FlexModul® ISO formats stands for modular flexibility and high quality of absolute tidy storage in cabinet modules

Trays and baskets can of course be equipped with limit switches for fixed withdrawals.



See special catalog:
Mobile benches and Carts
Trays and Baskets

Materials



ABS Grey



PC transparant



Wire basket, Coated

Compartment



Fixed dividers



Variable compartments



Removable dividers
for quick cleaning

Heights



50 mm



100 mm



200 mm

Størrelser



DIN-modul
600 x 400 mm



Half DIN-modul
300 x 400 mm

Cabinet shelves



Trays and baskets



Chemical and poison storage with wire baskets



Chemical and poison cabinets with pull-out

